


Trays With A Control Cutout



Features:

- **3/8" thick clear polycarbonate**
- Heavy-duty break-resistant
- Optional Top Drop Clamps attached
- Optional black PVC rim attached
- Recommended for chairs with fixed joysticks

Features:

- **3/8" thick melamine-impregnated fiberboard**
- Very durable melamine grey or black laminate surface
- Optional Top Drop Clamps attached
- Optional black PVC rim attached
- Recommended for chairs with fixed joysticks

CLEAR TRAYS WITH A CONTROL CUTOUT

	Clear Tray with Top Drop Clamps Attached				Clear Tray Only (no attaching hardware)		
	PVC Rim Right Side Cutout	PVC Rim Left Side Cutout	No Rim Right Side Cutout	No Rim Left Side Cutout	PVC Rim Right Side Cutout	PVC Rim Left Side Cutout	No Rim Cutout (Reversible)
21W x 18D	T13R12A	T13L12A	T13R02A	T13L02A	T13R10X	T13L10X	T13E00X
23 1/2 x 22	T15R12A	T15L12A	T15R02A	T15L02A	T15R10X	T15L10X	T15E00X
25 1/2 x 22	T16R12A	T16L12A	T16R02A	T16L02A	T16R10X	T16L10X	T16E00X
27 1/2 x 23 1/2	T17R12A	T17L12A	T17R02A	T17L02A	T17R10X	T17L10X	T17E00X

GREY TRAYS WITH A CONTROL CUTOUT

	Grey Tray with Top Drop Clamps Attached				Gray Tray Only (no attaching hardware)			
	PVC Rim Right Side Cutout	PVC Rim Left Side Cutout	No Rim Right Side Cutout	No Rim Left Side Cutout	PVC Rim Right Side Cutout	PVC Rim Left Side Cutout	No Rim Right Side Cutout	No Rim Left Side Cutout
21W x 18D	T23R12A	T23L12A	T23R02A	T23L02A	T23R10X	T23L10X	T23E00X	T23L00X
23 1/2 x 22	T25R12A	T25L12A	T25R02A	T25L02A	T25R10X	T25L10X	T25E00X	T25L00X
25 1/2 x 22	T26R12A	T26L12A	T26R02A	T26L02A	T26R10X	T26L10X	T26E00X	T26L00X
27 1/2 x 23 1/2	T27R12A	T27L12A	T27R02A	T27L02A	T27R10X	T27L10X	T27E00X	T27L00X

BLACK TRAYS WITH A CONTROL CUTOUT

	Black Tray with Top Drop Clamps Attached				Black Tray Only (no attaching hardware)			
	PVC Rim Right Side Cutout	PVC Rim Left Side Cutout	No Rim Right Side Cutout	No Rim Left Side Cutout	PVC Rim Right Side Cutout	PVC Rim Left Side Cutout	No Rim Right Side Cutout	No Rim Left Side Cutout
21W x 18D	T33R12A	T33L12A	T33R02A	T33L02A	T33R10X	T33L10X	T33E00X	T33L00X
23 1/2 x 22	T35R12A	T35L12A	T35R02A	T35L02A	T35R10X	T35L10X	T35E00X	T35L00X
25 1/2 x 22	T36R12A	T36L12A	T36R02A	T36L02A	T36R10X	T36L10X	T36E00X	T36L00X
27 1/2 x 23 1/2	T37R12A	T37L12A	T37R02A	T37L02A	T37R10X	T37L10X	T37E00X	T37L00X

Top Drop Clamps
position tray straight down from top.

Top Drop Clamps can be used on
Quantum TB3 arms by using our
TB3 Adapter Plates



Hook & Loop Straps


Pair, 2"x 24" 31091

Trays with Top Drop Clamps fit
a wheelchair with
Armpad-To-Armpad Width
(Outside Edges):

21W x 18D	from 16" to 18"
23 1/2W x 22D	from 18" to 20 1/2"
25 1/2W x 22D	from 20" to 22 1/2"
27 1/2W x 23 1/2D	from 22" to 24 1/2"



Pediatric



Trays with no Clamps fit a
wheelchair with
Seat Width
(Outside Edges):

21W x 18D	from 12" to 14"
23 1/2W x 22D	from 14" to 17"
25 1/2W x 22D	from 17" to 19"
27 1/2W x 23 1/2D	from 20" to 21"