TheraFit Premium Soft Goods

D Two-piece with sewn straps, hook and loop

center release with webbing pull loop

Hardware to attach to chair included

erafi

emium Soft Good

TheraFit Chest Strap lengths shown

here are taken from end to end, fully extended, without attaching

hardware. 1-Piece padded area

Medium 4" x 20", Large 4" x 24", X-Large 6" x 28", 2X-Large 8" x 32"

2-Piece padded area dimensions

are: Small 3" x 10" (each piece),

Medium 4" x 14", Large 4" x 18".

dimensions are: Small 3" x 16",

THERAFIT CHEST STRAPS

Features:

from skin irritation

C One-piece with sewn straps





